



Science.
Applied to Life.™

3M™ Lava™ Esthetic
Esthetic Fluorescent Full-Contour Zirconia

**Discover the next level
of full-contour esthetics.**

Natural esthetics. Now built in.



All you ever wanted in one crown: Ask your lab for 3M™ Lava™ Esthetic Zirconia – the high-quality cubic zirconia from 3M – and you will benefit from a strong and esthetic option for your crown and bridge cases.

Featuring a one-of-a-kind shading technology, Lava Esthetic zirconia combines the esthetics of glass ceramics with high full-contour strength.



For the first time, a zirconia that doesn't compromise on esthetics.



- Predictable, excellent shade match to VITA® Classic shades by built-in shade gradient.
- High translucency, optimized for esthetic full-contour restorations.
- First zirconia with inherent toothlike fluorescence for all shades.

Whether in daylight, sunshine or black light – restorations made of Lava Esthetic zirconia always look natural.





Sintered anterior crowns photographed on a fluorescent stump material.

Stronger than glass ceramics and other leading cubic zirconia.

With 800 MPa* strength and a high translucency optimized for esthetic full-contour crowns and bridges, Lava Esthetic zirconia solves typical glass ceramic cases with more confidence. With minimum 0.8 mm, it allows for reduced wall thickness versus IPS e.max® and with 12 mm² anterior (14 mm² posterior) also thinner connectors.

For cases with very limited space or parafunctional patients requiring ultimate strength, choose 3M™ Lava™ Plus Zirconia.

* 3-point bending strength according to ISO 6872:2015; qualified for Type II, class 4; indications: crowns, bridges with one pontic between two crowns, inlays, onlays and veneers.




Source: 3M Oral Care internal data. Data available upon request. See last page for contact information.

Indications:

- Crowns
- 3-unit bridges
- Inlays/onlays
- Veneers

Please refer to the instructions for use for details.

Antagonist friendly.




A 2016 study has shown that full-contour restorations made of Lava Esthetic zirconia are wear-friendly against opposing enamel.

Wear of enamel antagonist against incisal enamel and 3M™ Lava™ Esthetic Fluorescent Full-Contour Zirconia.

Source: Dr. J. O. Burgess, Dr. N. Lawson, University of Alabama at Birmingham School of Dentistry, U.S.A., Internal Report to 3M Oral Care, 2016.

Seating made easy.




Easy adjustment with fine-grain diamond bur and water cooling. Repolish adjusted area.

With your Lava Esthetic zirconia crown or bridge you receive a digital high-precision restoration from your lab. During try-in, handle crown gently and do not force it in place. If minor adjustments are required, Lava Esthetic zirconia can be more easily adjusted than conventional zirconia. Make sure to maintain a minimum wall thickness of 0.8 mm and re-polish adjusted areas.



For final cementation, 3M recommends using 3M™ RelyX™ 2 Unicem Automix Self-Adhesive Resin Cement offering high bond strength and ease of use without compromising reliability and esthetics. And, as there is no need to etch or prime, this easy-to-use cement will save you precious time and work steps.



Application of 3M™ RelyX™ Unicem 2 Automix Cement.

Please refer to the instructions for use for the proper cementation procedure.

3M™ RelyX™ Unicem 2 Automix Cement

Refills	Shade
56846	A2 Universal
56847	A3 Opaque
56848	Translucent

Each refill contains 1 syringe – 8.5 g,
15 mixing tips (regular), 5 mixing tips (wide), 5 endo tips

Ask your 3M sales representative for your copy of the

3M™ Lava™ Esthetic Fluorescent Full-Contour Zirconia Crown Quick Start Guide

or download from 3M.com/LavaEsthetic



3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000
USA
1-800-634-2249

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-888-363-3685

3M.com/LavaEsthetic

3M, 3M Science. Applied to Life., Lava and RelyX are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada. All other trademarks are owned by other companies. Please recycle. Printed in U.S.A. © 3M 2016. All rights reserved. 70-2013-0831-2